


DIAGNOSTIC WEB VIEWER

iQ-X is a powerful iQ-WEB plug-in for browsing, reading, and processing DICOM images using Internet Explorer®.

iQ-X provides an interactive diagnostic web viewer with all clinically required image processing features such as measurements, stack mode, cine mode, and study comparisons.

iQ-X allows additional series to be cached in the background while images are being viewed. The diagnostic image processing tools available are based on the leading DICOM viewer iQ-VIEW.

Web-Based Reporting: With the optional iQ-X reporting module, physicians can create DICOM Structured Reports online from anywhere in the world.



iQ-X WEB VIEWER

- View any kind of images, Structured Reports, and Encapsulated PDFs**
- Thumbnail previews for easy selection of series via drag & drop
- Displays scoutlines of cross-sectional images
- Center/window, color remappings
- Stack mode/cine mode
- Measurements (distance, ROI, angle, modify/delete options)
- Basic annotation of images
- Support of look-up tables, bitmap overlays, DICOM shutter sequences
- Compare multiple studies
- Easy synchronization of series, e.g. for comparing CT scans
- Web-based creation of DICOM Structured Reports with Latin-1 support (optional)
- Background caching
- JPEG / JPEG LS / JPEG 2000 support (lossless/lossy)

LICENSING

- Unlimited number of iQ-X clients
- iQ-X is available with or without report editor

CERTIFICATION

- CE 0482 and FDA 510(k)***



* Detailed information on iQ-X's system requirements is available in the iQ-X Administration Guide or the iQ-SYSTEM PACS Hardware Purchasing Guide.

** PDF reader is required.

*** iQ-X is a component of iQ-SYSTEM PACS 1.4 (FDA UDI B403IQS140).

